



Application of Aerogel for Prosthetic Liners

Carolina I. Ragolta
Mentors: Luke Roberson & Michele Birmele
NASA MUST Intern, Kennedy Space Center
August 11, 2011



Introduction





About Me

- From Miami, FL
- Junior at Harvard College
 - Major: BME
 - Minor: Economics
- 2nd KSC internship through the NASA MUST Program







Summer Experiences






Background



Prosthetic Liners

- Interface between residual limb and prosthetic
- Provide comfort and control
- Many different materials
 - Cotton
 - Wool
 - Synthetics
 - Silicone or urethane gel

Existing Liner Problems


- Tradeoff between comfort and durability
- Very hot
- Airtight – sweat pools
 - Maintenance
 - Itching
 - Odor
 - Discomfort
 - Friction (leads to injury)



KSC
ENGINEERING


Aerogel Properties

- Hydrophobic
- Breathable
- Flexible
- Environmentally friendly
- Non-toxic
- Used in hot- and cold-weather clothing
- Blanket used for project encased in nylon



KSC
ENGINEERING

How Aerogel is Made



KSC
ENGINEERING

Images from www.aerogel.org

Design Requirements

- Perspiration control
- Stay cool in FL heat
- Distribute load evenly
- Prevent injury
- Prevent skin irritation
- Accommodate volume fluctuation in residual limb
- Easy to don and doff



KSC
ENGINEERING

Test Methods

KSC
ENGINEERING

Moisture Vapor Permeability


- No ISO standards for interface materials
- Literature search: one study using distilled water
- Adapted test for wound dressings
 - BS EN 13726-1:2002
 - Artificial sweat solution from ISO 3160-2:2003



KSC
ENGINEERING

Moisture Vapor Permeability


- Samples prepared in the Prototype Lab
- Circular samples to fit test assembly (4 cm diameter)
- Encased in nylon to prevent skin irritation
- Sealed with waterproof First Aid tape



KSC
ENGINEERING

Moisture Vapor Permeability


- Test assembly with sample and fluid weighed
- Incubated inverted
 - 24 hours
 - 37°C, 5% relative humidity
- Test assembly removed and reweighed
 - Difference in mass is the fluid that has transpired



KSC
ENGINEERING

Moisture Vapor Permeability


- Aerogel beads
 - Difficult to use
 - Concerns about force transmission
- Comparison of Pyrogel® (2.0 mm and 6.0 mm) and 10.0 mm Spaceloft®
 - Spaceloft® difficult to use
 - Testing continued with the two Pyrogel® thicknesses



KSC
ENGINEERING

Compression


- 1.5" square samples of Pyrogel®
- Tested on Instron
 - Max load: 1000 lbf
 - Compression rate: 0.1" per minute
- Thickness measured before and after



KSC
ENGINEERING

Biofilm Formation

- 12 mm diameter coupons of Pyrogel®
- Tested on CDC Biofilm Reactor
 - ASTM E 2562-07
 - Challenge organism: *Pseudomonas aeruginosa*
 - 24 hour batch mode
 - 24 hour CSTR mode



KSC
ENGINEERING

Results and Discussion


KSC
ENGINEERING

Moisture Vapor Permeability

Material	Moisture Vapor Lost (g)	Moisture Vapor Lost (%)
Nanogel®	-12.8	84.7
	-12.3	86.9
	-12.9	85.7
2.0 mm Pyrogel®	-13.4	88.4
	-5.19	34.8
	-6.55	43.8
6.0 mm Pyrogel®	-2.43	16.1
	-3.48	23.0
	-10.1	67.3
10.0 mm Spaceloft®	-6.27	41.4

KSC
ENGINEERING

Moisture Vapor Permeability



KSC
ENGINEERING


Why such high variance?

- Non-uniform performance of the test fixtures
- Leaking was noted in two trials
 - Last two trials
 - Where failure was noted, data points were eliminated
- May have leaked in other trials
 - Fluid may have evaporated before it was detected
 - Did not keep track of which test fixture went with which sample for each trial

KSC
ENGINEERING

Comparison with Existing Materials

- Existing materials allow little to no moisture to pass through
- Further testing is needed
 - No test standards
 - Only study published has used distilled water
 - Comparison with GORE-TEX®, also breathable



KSC
ENGINEERING


Chart from Hachisuka, K., Matsushima, Y., Ohmine, S., Shitama, H., and Shinkoda, K., "Moisture permeability of the total surface bearing prosthetic socket with a silicone liner: is it superior to the patella-tendon bearing prosthetic socket?" Journal of UOEH, Vol. 23, No. 3, 2001, pp. 225-232.

Compression

Material	Initial Thickness (mm)	Compression at 1000 N (mm)	Strain at 1000 N (%)	Loss in Thickness (mm)	Loss in Thickness (%)
2.0 mm Pyrogel®	2.11	1.74	77.4	-0.432	-20.5
	2.18	1.32	51.9	-0.483	-22.1
	2.39	1.45	44.6	-0.483	-20.2
6.0 mm Pyrogel®	5.94	3.78	64.5	-0.356	-5.98
	6.02	3.64	60.2	-0.737	-12.2
	6.22	3.47	52.2	-0.889	-14.3

KSC
ENGINEERING


Compression



KSC
ENGINEERING

Comparison with Existing Materials

- Gait analysis: maximum axial force is between 800 and 900 N
- Previous studies: 550 N
- This study: over 4400 N
- Performed comparable to or better than existing materials, even under eight times the load




KSC
ENGINEERING

Conclusions

KSC
ENGINEERING

Suitability for Prosthetic Liners

- Breathability
 - Permeable to vapor
 - Further testing needed with existing materials
- Load bearing
 - Performed comparably under higher loads
 - Further testing needed with existing materials



KSC
ENGINEERING

Future Development

- Biofilm testing
- Comparisons with existing materials
- Uniformly performing test fixtures
- Repeated or cyclic load bearing tests
- Friction load bearing tests
- Liner prototype



KSC
ENGINEERING

Acknowledgements

- | | |
|-------------------|----------------|
| • Luke Roberson | • Jeff Sampson |
| • Tim Griffin | • Dean Lewis |
| • Dionne Jackson | • Brian Taylor |
| • Annie Caraccio | • Don Doerr |
| • Dan Woodard | • Barry Slack |
| • Wes Johnson | • Ron Woods |
| • Michele Birmele | • NE-L |
| • Megan Morford | |



THANK YOU SO MUCH!

KSC
ENGINEERING

Questions

KSC
ENGINEERING